



**better analysis counts**

## **Фокусирующая поликапиллярная оптика Техническое описание**

Архангельск (8182)63-90-72  
Астана +7(7172)727-132  
Белгород (4722)40-23-64  
Брянск (4832)59-03-52  
Владивосток (423)249-28-31  
Волгоград (844)278-03-48  
Вологда (8172)26-41-59  
Воронеж (473)204-51-73  
Екатеринбург (343)384-55-89  
Иваново (4932)77-34-06  
Ижевск (3412)26-03-58  
Казань (843)206-01-48

Калининград (4012)72-03-81  
Калуга (4842)92-23-67  
Кемерово (3842)65-04-62  
Киров (8332)68-02-04  
Краснодар (861)203-40-90  
Красноярск (391)204-63-61  
Курск (4712)77-13-04  
Липецк (4742)52-20-81  
Магнитогорск (3519)55-03-13  
Москва (495)268-04-70  
Мурманск (8152)59-64-93  
Набережные Челны (8552)20-53-41

Нижний Новгород (831)429-08-12  
Новокузнецк (3843)20-46-81  
Новосибирск (383)227-86-73  
Орел (4862)44-53-42  
Оренбург (3532)37-68-04  
Пенза (8412)22-31-16  
Пермь (342)205-81-47  
Ростов-на-Дону (863)308-18-15  
Рязань (4912)46-61-64  
Самара (846)206-03-16  
Санкт-Петербург (812)309-46-40  
Саратов (845)249-38-78

Смоленск (4812)29-41-54  
Сочи (862)225-72-31  
Ставрополь (8652)20-65-13  
Тверь (4822)63-31-35  
Томск (3822)98-41-53  
Тула (4872)74-02-29  
Тюмень (3452)66-21-18  
Ульяновск (8422)24-23-59  
Уфа (347)229-48-12  
Челябинск (351)202-03-61  
Череповец (8202)49-02-64  
Ярославль (4852)69-52-93



**better analysis counts**

# Polycapillary Optics

**FOR MICRO X-RAY FLUORESCENCE  
AND X-RAY DIFFRACTION**



## **Improved Performance with Polycapillary Optics:**

State-of-the-art optic captures a large solid angle from an X-ray source and redirects the beam into a micron-sized focal spot or highly collimated beam.

### **Features:**

- Orders of magnitude flux gain from micron-sized spot
- Integration with compact, low-power sources provides flux equivalent to rotating anode sources
- Broad spectral bandwidth: 50 eV- 50 KeV
- Point-to-point focusing beam
- Point-to-parallel beam
- Custom designed enclosure available

### **Benefits:**

- Extremely high flux density
- Increased spatial resolution
- Ideal for analysis of irregular shaped, unprepared, or low-Z samples

### **Advantages over Electron Probe X-ray Analysis:**

- Enhanced detection sensitivity
- No special sample preparation needed
- Operation in air

### **Applications:**

- Evaluation of small features
- Elemental mapping
- Film/plating thickness
- Detection of micro-contamination
- Multi-layered coatings for advanced circuit boards
- Small particle analysis
- Forensics
- Powder XRD
- Single crystal XRD

# XRF and XRD with Polycapillary Optics

**Custom Solutions.** In XRF applications, focusing polycapillary optics deliver micron-sized focal spots and beam intensity superior to pinhole based micro XRF. Several orders of magnitude flux gain are possible, depending on the experimental geometry. In diffraction applications, polycapillary optics offer a highly collimated beam. Below are typical examples of polycapillary optics. Customized polycapillary optics based on customer requirements are available to meet specific instrument requirements.

## Focusing Optics ( $\mu$ XRF)

|                                                                               |      |      |      |      |     |     |     |
|-------------------------------------------------------------------------------|------|------|------|------|-----|-----|-----|
| OFD (mm)                                                                      | 2    | 4    | 9    | 20   | 50  | 100 | 200 |
| Focal spot size ( $\mu$ m, FWHM, Mo Ka)                                       | 8    | 15   | 25   | 45   | 100 | 180 | 300 |
| Intensity gain (vs a pinhole collimation of same size, 100mm from the source) | 6000 | 3000 | 2200 | 1200 | 400 | 160 | 70  |

### Note:

1. The above performance parameters relate to 50  $\mu$ m X-ray source.
2. The IFD of the optic is 20mm.

## Collimating/Parallel-beam Optics (XRD/WDS)

|                           |     |    |     |     |     |     |     |     |
|---------------------------|-----|----|-----|-----|-----|-----|-----|-----|
| Output beam diameter (mm) | 0.5 | 1  | 2   | 3   | 4   | 6   | 10  | 15  |
| Intensity gain            | 12  | 45 | 130 | 250 | 370 | 470 | 680 | 850 |

### Note:


1. The above performance parameters are for a 50 $\mu$ m X-ray source at 8keV.
2. The IFD of the optic is 18mm and the output divergent angle is 0.2 degree.
3. Intensity gain determined by comparison to pinhole configuration at same distance from source and has the same beam size and divergent angle.

## Slightly-focusing Optics (XRD/XRF)


|                                  |      |     |     |     |
|----------------------------------|------|-----|-----|-----|
| Output convergent angle (degree) | 0.25 | 0.5 | 1   | 2   |
| Intensity gain                   | 25   | 80  | 260 | 820 |

### Note:

\* 50 $\mu$ m X-ray source IFD: 18mm OFD: 140mm Energy: 8keV Focal spot size: 0.5mm, FWHM.




**Figure 3** Peak shift as function of sample displacement for parallel-beam geometry and Bragg-Brentano geometry.



**Figure 1**

Comparison of  $\mu$ XRF spectra generated using a focusing polycapillary optic and a pinhole aperture. The spectra shown is of an air particulate sample  $\sim 50\mu$ m in diameter. Mo excitation (40kV, 20W).



**Figure 2**

Spectrum of NIST SRM 1832 standard XRF sample using a polycapillary optic. Mo excitation (40kV, 12W).

**Table 1**

Picogram detection limits are possible using the polycapillary optic for this thin film SRM sample.

| Element                 | K   | Ti  | Fe   | Zn   | Pb   |
|-------------------------|-----|-----|------|------|------|
| Limit of Detection (pg) | 4.1 | 1.5 | 0.57 | 0.28 | 0.52 |



**better analysis counts**



**better analysis counts**

Архангельск (8182)63-90-72  
Астана +7(7172)727-132  
Белгород (4722)40-23-64  
Брянск (4832)59-03-52  
Владивосток (423)249-28-31  
Волгоград (844)278-03-48  
Вологда (8172)26-41-59  
Воронеж (473)204-51-73  
Екатеринбург (343)384-55-89  
Иваново (4932)77-34-06  
Ижевск (3412)26-03-58  
Казань (843)206-01-48

Калининград (4012)72-03-81  
Калуга (4842)92-23-67  
Кемерово (3842)65-04-62  
Киров (8332)68-02-04  
Краснодар (861)203-40-90  
Красноярск (391)204-63-61  
Курск (4712)77-13-04  
Липецк (4742)52-20-81  
Магнитогорск (3519)55-03-13  
Москва (495)268-04-70  
Мурманск (8152)59-64-93  
Набережные Челны (8552)20-53-41

Нижний Новгород (831)429-08-12  
Новокузнецк (3843)20-46-81  
Новосибирск (383)227-86-73  
Орел (4862)44-53-42  
Оренбург (3532)37-68-04  
Пенза (8412)22-31-16  
Пермь (342)205-81-47  
Ростов-на-Дону (863)308-18-15  
Рязань (4912)46-61-64  
Самара (846)206-03-16  
Санкт-Петербург (812)309-46-40  
Саратов (845)249-38-78

Смоленск (4812)29-41-54  
Сочи (862)225-72-31  
Ставрополь (8652)20-65-13  
Тверь (4822)63-31-35  
Томск (3822)98-41-53  
Тула (4872)74-02-29  
Тюмень (3452)66-21-18  
Ульяновск (8422)24-23-59  
Уфа (347)229-48-12  
Челябинск (351)202-03-61  
Череповец (8202)49-02-64  
Ярославль (4852)69-52-93